

POSTER PRESENTATION

Cine interleaved sequences enabled imaging of mice on clinical 3T MRI and analysis of their cardiac function after myocardial infarction

Alexandre Belin^{1,2*}, Vincent Braunersreuther^{1,2}, Fabrizio Montecucco^{1,2}, Benedicte M Delattre^{1,2}, Francois Mach^{1,2}, Jean Vallee^{1,2}

From 15th Annual SCMR Scientific Sessions Orlando, FL, USA. 2-5 February 2012

Background

With the poor availability of small animal dedicated MRI, it is of great interest and challenge to use clinical MRI to image and analyze the cardiac function of mice. This would of course be advantageous for translational research and also enable the use of up-to-date sequences already implemented in clinical routine. The aim is to study the time evolution of cardiac function in mice with a myocardial infarct on a clinical MRI.

Methods

C57BL/6 mice (n=4) were submitted in vivo to left coronary artery permanent ligature and compared and with compared with non-operated mice (n=4). The mice were imaged on a clinical Siemens 3T MRI using an "interleaved" sequence constructed from an ECG-triggered turboflash cine sequence, combining two acquisitions shifted in time yielding an effective time resolution of 6.8 ms and 20-26 phases per heart beat with following parameters: FOV 111 mm, in-plane resolution 257 µm, slice thickness 1 mm, TE/TR 6.2/13.5 ms, flip angle 30°. A soft-thresholding of the temporal Fourier coefficients was used to further denoise the images. The mice were scanned at 24h and 22 days after coronary ligation.

Results

High quality images suitable for endocardial contouring were obtained for all the animals. The systolic and diastolic volumes, as well as the ejection fraction, are reported in Table 1. The anterior wall of all the operated animals was akenetic as shown on Figure 1. The

¹Radiology, HUG, Geneve 14, Switzerland

Full list of author information is available at the end of the article



© 2012 Belin et al; licensee BioMed Central Ltd. This is an open access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/2.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

ejection fraction was significantly decreased in the infarcted group at 24h and 22 days by comparison to the control group (p=0.006 and p=0.003). There was also a trend for the global function to decrease from 24h to 22 days.

Table 1 Cardiac properties of the mice

	ESV (ul)	EDV (ul)	EF
INF 24h			
EVAMR135	65.26	88.12	0.26
EVAMR90	51.27	69.73	0.26
EVAMR174	43.66	70.13	0.38
EVAMR175	51.27	81.90	0.37
MEAN	52.87	77.47	0.32
SD	9.01	9.07	0.07
INF 22j			
EVAMR135	289.57	356.57	0.19
EVAMR90	208.26	241.62	0.14
EVAMR174	61.75	91.41	0.32
EVAMR175	66.65	109.09	0.39
MEAN	156.56	199.67	0.26
SD	111.71	124.23	0.12
NON INF			
P1	21.13	56.86	0.63
P2	33.92	69.12	0.51
P3	36.83	79.46	0.54
P4	47.59	90.64	0.47
MEAN	34.87	74.02	0.54
SD	10.88	14.43	0.07

Open Access



Conclusions

We demonstrated a robust protocol to study cardiac functions in mice with a myocardial infarct using a clinical 3T MRI. This methodology has a strong potential to study the effect of treatment in rodents.

Funding

This work was partially supported by the Swiss National Science Foundation (grant PP00P2-123438).

Author details

¹Radiology, HUG, Geneve 14, Switzerland. ²Faculty of medicine, University of Geneva, Geneva, Switzerland.

Published: 1 February 2012

doi:10.1186/1532-429X-14-S1-P60 **Cite this article as:** Belin *et al.*: Cine interleaved sequences enabled imaging of mice on clinical 3T MRI and analysis of their cardiac function after myocardial infarction. *Journal of Cardiovascular Magnetic Resonance* 2012 14(Suppl 1):P60.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

BioMed Central

Submit your manuscript at www.biomedcentral.com/submit