

ORAL PRESENTATION



Right ventricular kinetic energy: 4D flow MRI analysis of healthy volunteers and repaired Tetralogy of Fallot

Daniel Jeong^{1*}, Alejandro Roldán-Alzate², Christopher J Francois¹

From 17th Annual SCMR Scientific Sessions New Orleans, LA, USA. 16-19 January 2014

Background

Repaired Tetralogy of Fallot (rTOF) is currently monitored with Cardiovascular MR (CMR) to measure right ventricle (RV) volumes and function. However, these indices do not always reflect performance status or predict success following repair of pulmonary insufficiency. RV kinetic energy (KERV) may provide a more robust measure of RV efficiency than RV end diastolic volume (EDV) index and function. The aim of this study was to demonstrate the feasibility of calculating KERV with 4D Flow MRI and assess differences between rTOF and healthy volunteers.


Methods

Ten subjects with rTOF and 9 healthy volunteers were scanned according to an IRB approved protocol. 4D flow MRI (Phase Contrast Vastly Undersampled Isotropic Projection Reconstruction PCVIPR) was performed on clinical MRI scanners (GE Healthcare, Waukesha, WI) after administration of gadolinium based contrast agents. PC VIPR parameters: FOV = 260-320 mm3, isotropic 1.3 mm spatial resolution, TR/TE = 8.8-10.9/2.8-3.7 ms (1.5T), and 6.2-3.7/2.0-2.2 ms (3.0T), Venc = 40-400 cm/s, scan time: 9-17 minutes using respiratory and retrospective ECG gating. Post processing was done using Mimics (Materialise, Leuven, Belgium) and MatLab (The MathWorks, Natick, MA) to measure KERV and main pulmonary artery flow (QMPA) throughout the cardiac cycle. The KE of a voxel of blood was calculated from 1/2 mv^2, where the mass (m) represents the voxel volume multiplied by the density of blood (1.025 g/ml) and the velocity (v) of each

¹Department of Radiology, University of Wisconsin, Madison, Wisconsin, USA Full list of author information is available at the end of the article voxel was determined from 4D flow MRI. KERV was determined from the sum of the KE of the voxels within the segmented RV volume at each phase of the cardiac cycle. Differences in peak systolic and diastolic KERV and the QMPA/KERV ratio, an indicator of RV-MPA coupling, between groups were assessed using the Wilcoxon rank-sum test.


Results

Three distinct KE peaks were observed in all subjects. Peak systolic (p = 0.0076) and diastolic (p = 0.0003 peak 1, p = 0.013 peak 2) were higher in rTOF than in healthy volunteers (Figures 1,2). The QMPA/KERV ratio was lower in rTOF (0.92 +/- 0.68 ml/cycle*mJ) than in healthy volunteers (6.92 +/- 4.16, p = 0.0002).





© 2014 Jeong et al.; licensee BioMed Central Ltd. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/2.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited. The Creative Commons Public Domain Dedication waiver (http:// creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated.



Conclusions

Time-resolved KERV was measured in rTOF and healthy volunteers using 4D flow MRI. KERV was higher and QMPA/KERV ratio was lower in rTOF than in healthy volunteers, indicative of greater inefficiency in RV function to generate the same cardiac output. Future studies are needed to determine if changes in RV KE provide earlier evidence of RV dysfunction and need for reintervention than standard CMR measurements.

Funding

NIH grant R01HL07226. GE Healthcare for research support.

Authors' details

¹Department of Radiology, University of Wisconsin, Madison, Wisconsin, USA. ²Department of Medical Physics, University of Wisconsin, Madison, Wisconsin, USA.

Published: 16 January 2014

doi:10.1186/1532-429X-16-S1-O46

Cite this article as: Jeong *et al.*: **Right ventricular kinetic energy: 4D flow MRI analysis of healthy volunteers and repaired Tetralogy of Fallot.** *Journal of Cardiovascular Magnetic Resonance* 2014 **16**(Suppl 1):O46.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

BioMed Central

Submit your manuscript at www.biomedcentral.com/submit