

WALKING POSTER PRESENTATION

Open Access

Fast, whole-heart, free-breathing 3D T₂ mapping at 3T with application to myocardial edema imaging

Hsin-Jung Yang^{1*}, Jianing Pang¹, Behzad Sharif¹, Avinash Kali¹, Xiaoming Bi², Ivan Cokic¹, Debiao Li¹, Rohan Dharmakumar¹

From 18th Annual SCMR Scientific Sessions
Nice, France. 4-7 February 2015

Background

Cardiac MRI (CMR) T₂ mapping is a proven method for myocardial edema detection. However, the current approach requires multiple breath-holds and can take nearly 10-15 minutes to complete. Faster acquisitions could potentially improve patient comfort and cost-effectiveness of CMR exams. The objective of this study is to develop and test a free-breathing, three-dimensional, cardiac MR approach which can yield fast, accurate, T₂ maps of the whole left ventricle at 3T.

Methods

We developed an ECG-triggered, free-breathing, T₂-prepared, three-dimensional gradient-echo acquisitions with different echo times (0, 24, 55 ms) with near perfect navigator efficiency on a clinical 3T system. The proposed approach was tested and validated in ex-vivo porcine hearts, healthy volunteers and canines with reperfused acute myocardial infarction (rAMI). On the basis of the navigator signals, images were corrected for respiratory motion and were fit to a mono-exponential function to derive T₂ maps of the whole left-ventricular myocardium.

Results

Ex-vivo myocardial T₂ values of the proposed approach (3D FB MoCo) were not different from standard 2D approaches (all $p < 0.05$): 48.7 ± 0.9 ms (3D FB MoCo) vs. 48.2 ± 0.6 ms (2D spin echo) and 47.5 ± 0.8 ms (2D T₂-prepared bSSFP (T₂-prep bSSFP)). In healthy volunteers, compared to 3D FB MoCo and 2D BH, myocardial T₂

maps, 3D FB Non-MoCo T₂ myocardial maps showed longer T₂ values ($p < 0.05$), larger coefficient-of-variations (COV) in T₂ ($p < 0.05$), and lower image quality ($p < 0.05$). Conversely, the mean and COV in myocardial T₂ and image quality of 2D BH and 3D FB MoCo T₂ were not different ($p = 0.99$, $p = 0.74$, $p = 0.14$, respectively). In canines with rAMI, edema volumes measured from 2D BH and 3D FB MoCo T₂ maps were closely correlated (both $R^2 = 0.97$ and $p < 0.05$). In Bland-Altman analysis, mean T₂ of edematous and remote zones and edema volumes were within the limits of agreement (bias in T₂ = 0.4 ms and edema volume = 0.9%).

Conclusions

The proposed free-breathing, three-dimensional T₂ mapping approach at 3T enabled whole-heart acquisitions within 5 minutes with an accuracy in T₂ not different from that of the state-of-the-art breath-held T₂ mapping approach.


Funding

This work was supported in part by National Heart, Lung, and Blood Institute HL091989 (R.D.). The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Heart, Lung, and Blood Institute or the National Institutes of Health.

Authors' details

¹Cedars Sinai Medical Center, Los Angeles, CA, USA. ²Siemens Healthcare, Los Angeles, CA, USA.

¹Cedars Sinai Medical Center, Los Angeles, CA, USA
Full list of author information is available at the end of the article



Published: 3 February 2015

doi:10.1186/1532-429X-17-S1-Q125

Cite this article as: Yang *et al.*: Fast, whole-heart, free-breathing 3D T2 mapping at 3T with application to myocardial edema imaging. *Journal of Cardiovascular Magnetic Resonance* 2015 **17**(Suppl 1):Q125.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at
www.biomedcentral.com/submit

