WORKSHOP PRESENTATION

Open Access

Four-dimensional velocity encoded MRI improves blood flow quantification in patients with semilunar valve stenosis

Sarah Nordmeyer*, Eugenie Riesenkampff, Daniel Messroghli, Felix Berger, Titus Kuehne

From 15th Annual SCMR Scientific Sessions Orlando, FL, USA. 2-5 February 2012

Background

Stenoses might cause complex flow patterns, which are sometimes difficult to assess quantitatively with standard two-dimensional (2D) VEC MRI. We sought to evaluate the use of four-dimensional (4D) velocity encoded magnetic resonance imaging (VEC MRI) for blood flow quantification in patients with semilunar valve stenosis.

Methods

Peak velocities (Vmax) and stroke volumes (SV) were quantified by 2D and 4D VEC MRI in volunteers (n=7) and patients with semilunar valve stenosis (n=18). Measurements were performed immediately above the aortic or pulmonary valve with both techniques (=level 1) and, additionally, at further predefined planes in the ascending aorta and in the pulmonary trunk within the 4D dataset. 4D VEC MRI streamline analysis identified "individualized" planes of highest flow velocity (4Dmaxtargeted) for further measurements. In patients, Vmax was also measured by Doppler-echocardiography.

Results

In patients, 4D VEC MRI (2.7m/s) showed comparable Vmax to Doppler-echocardiography (2.8m/s) and significantly higher Vmax than 2D VEC MRI (2.4m/s) (p<0.03) at level 1. 4Dmax-targeted revealed highest Vmax values (3.1m/s). Correlations of MR-derived peak velocities with Doppler-echocardiography were r=0.62 for 2D, r=0.67 for 4D at level 1 and r=0.80 for 4Dmax-targeted. 4D showed higher agreement with Doppler-echocardiography than did 2D VEC MRI at level 1 (p=0.039).

SV at level 1 were comparable between both techniques. SV measurements at different anatomical levels in the ascending aorta showed a significantly larger variance in patients with complex flow patterns than in volunteers with laminar flow patterns (p=0.004).

Conclusions

4D VEC MRI improves MR-derived quantification of peak flow velocities in patients with semilunar valve stenosis. Additionally, 4D VEC MRI demonstrates the potential for SV underestimation in complex flow.

Published: 1 February 2012

doi:10.1186/1532-429X-14-S1-W10

Cite this article as: Nordmeyer *et al*: Four-dimensional velocity encoded MRI improves blood flow quantification in patients with semilunar valve stenosis. *Journal of Cardiovascular Magnetic Resonance* 2012 14(Suppl 1): W10.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at www.biomedcentral.com/submit



Department of Congenital Heart Disease/Pediatric Cardiology, Deutsches Herzzentrum Berlin, Berlin, Germany



Table 1 Comparison between 2D (level 1) and 4D (level 1 and 4D max-targeted) derived peak velocity (+Doppler-echocardiography) in the ascending aorta in patients with aortic valve stenosis and in the pulmonary trunk in patients with pulmonary valve stenosis

| Patients | Vmax (m/s) ECHO | Vmax (m/s) 2D Level 1 | Vmax (m/s) 4D Level 1 | Vmax (m/s) 4D max-targeted |
|-----------|-----------------|-----------------------|-----------------------|----------------------------|
| Mean | 2.78 | 2.40 | 2.74 | 3.12 |
| SD | 0.68 | 0.54 | 0.64 | 0.63 |
| p-value* | 0.009 | - | 0.025 | < 0.001 |
| p-value** | - | <0.01 | >0.05 | <0.05 |

Vmax = peak velocity, 2D = two-dimensional, 4D = four-dimensional, SD = standard deviation, Echo = Doppler-echocardiography, * = compared to 2D, ** = compared to Echo




Figure 1 Assessment of peak velocity using flow visualization by 4D VEC MRI Figure: Velocity information in the left ventricular outflow tract of a patient with bicuspid aortic valve. A) Coronal axis of the left ventricular outflow tract, aortic valve and ascending aorta with in-plane visualization of velocity information. Dark red (black arrow) visualizes site of peak velocity. B) 4Dmax-targeted plane for blood flow quantification with colour coded flow direction.