

ORAL PRESENTATION

Open Access

Towards reliable myocardial blood-oxygen-level-dependent (BOLD) CMR using late effects of regadenoson with simultaneous ^{13}N -ammonia PET validation in a whole-body hybrid PET/MR system

Hsin-Jung Yang^{1*}, Damini Dey¹, Jane M Sykes², John Butler², Behzad Sharif¹, Debiao Li¹, Sotirios Tsaftaris³, Piotr Slomka¹, Frank S Prato², Rohan Dharmakumar¹

From 19th Annual SCMR Scientific Sessions
Los Angeles, CA, USA. 27-30 January 2016

Background

BOLD CMR is a non-contrast approach for examining myocardial perfusion but despite major technical advancements to date, its reliability remains weak. A key reason for this is the unpredictable cardiac motion during stress, which can lead to pronounced artifacts that confound/mask the true BOLD signal changes during hyperemia. Recently, regadenoson has become the vasodilator of choice owing to greater patient tolerability and ease of use. We hypothesized that at 10-mins post regadenoson administration (p.r.a), (a) BOLD CMR artifacts at stress are markedly reduced compared to those conventionally acquired at 2-mins p.r.a; and (b) that myocardial perfusion reserve (MPR) remains greater than 2.0 and is highly correlated with the BOLD effects estimated from T_2 maps.

Methods

Canines ($n = 7$) were studied in a PET/MR system. MR acquisitions were used to generate short-axis 2D T_2 maps; and the PET acquisitions following ^{13}N -ammonia infusion were used to quantify myocardial blood flow (MBF). Initially, 2D T_2 maps and PET signals were acquired at rest. Subsequently, regadenoson ($2.5 \mu\text{g}/\text{kg}$) was administered. T_2 maps were acquired at 2- and 10-mins p.r.a and PET signals were acquired at 10-mins p.r.a. Standard deviation (s) of myocardial T_2 values was measured at rest, 2- and 10-mins p.r.a from T_2 maps and were used to determine *Myocardial BOLD Variability* (MBV, defined as $sT_2(\text{stress})/sT_2(\text{rest})$) at 2- and 10-min p.r.a. Similarly, using

the mean T_2 values, *Myocardial BOLD Response* (MBR, defined as $T_2(\text{stress})/T_2(\text{rest})$) was computed at 10-mins p.r.a. PET images were analyzed with qPET software to determine MBF and MPR at rest and 10-mins p.r.a and were regressed against MBR

Results

A box-plot of observed MBV at 2- and 10-mins p.r.a (and at rest, for reference), along with representative T_2 maps are shown in Fig. 1. Note the extensive artifacts present in the T_2 map at 2 min, which are absent in the T_2 maps acquired at rest and 10-mins p.r.a. MBV was significantly larger at 2-mins p.r.a (1.6 ± 0.9) compared to 10-mins p.r.a (1.0 ± 0.3) and rest (1.0); $p < 0.05$ for both. Representative MBF at rest and 10-mins p.r.a are shown in Fig. 2A. MBF at 10-min p.r.a ($1.8 \pm 0.9 \text{ ml}/\text{g}/\text{min}$) was significantly higher than at rest ($0.6 \pm 0.3 \text{ ml}/\text{g}/\text{min}$), $p < 0.05$ (Fig. 2B). Mean MPR at 10-min p.r.a was 3.0. Corresponding BOLD images (T_2 maps) are shown in Fig. 2C. Myocardial T_2 at 10-min p.r.a ($40.4 \pm 1.7 \text{ ms}$) was significantly higher than at rest ($37.1 \pm 2.0 \text{ ms}$), $p < 0.05$ (Fig. 2D). MBR was strongly correlated with MPR ($R = 0.7$, $p < 0.05$, Fig. 2E)

Conclusions

Myocardial BOLD images acquired at 10-min p.r.a (compared to 2-min p.r.a) can be free of image artifacts. MPR at 10-mins p.r.a can be consistently higher than 2.0 and is strongly correlated with MBR. These data support that delayed acquisition of BOLD CMR post regadenoson administration is a viable means for increasing the reliability of cardiac BOLD. The clinical utility of this approach remains to be evaluated in human subjects.

¹Cedars Sinai Medical Center, Los Angeles, CA, USA
Full list of author information is available at the end of the article




Figure 1 Myocardial BOLD Variability at Stress Relative to Rest. Box-plot of Myocardial BOLD Variability (σT_2 (stress) / σT_2 (rest)) computed from T_2 maps acquired at 2 min and 10 min post regadenoson administration (p.r.a) is shown. The mean and range of Myocardial BOLD Variability are shown under various conditions. Large Myocardial BOLD Variability is observed at 2 mins p.r.a. compared to rest and is markedly reduced at 10 mins p.r.a. Representative images acquired at the various conditions are shown for reference. Note the marked T_2 inhomogeneity in T_2 image at 2 min p.r.a. compared to rest and 10 min p.r.a. σT_2 denotes the standard deviation of myocardial T_2 .




Figure 2 ^{13}N -Ammonia PET Myocardial Blood Flow and BOLD Response at 10 mins post regadenoson administration (p.r.a.). Representative rest and stress (10 mins p.r.a) short-axis PET Myocardial Blood Flow (MBF) and myocardial BOLD T_2 pages are shown in panels A and C. Both PET and BOLD images shown significant increase in MBF and BOLD response, respectively, at 10 mins p.r.a. compared to rest. Box plot of rest and stress MBF and myocardial T_2 across all animals are shown in panels B and D, respectively. Mean increase in MBF by a factor of 3.0 and a 9% T_2 elevation were observed on the PET and T_2 maps acquired at 20 min p.r.a. relative to rest. Results from regression analysis showed good correlation between PET myocardial perfusion reserve (MBF(stress)/MBF(rest)) and Myocardial BOLD Response ($R = 0.7$, $p < 0.05$; panel E).

Authors' details

¹Cedars Sinai Medical Center, Los Angeles, CA, USA. ²Lawson Health Research Institute, London, ON, Canada. ³Institute for Advanced Studies Lucca, Lucca, Italy.

Published: 27 January 2016

doi:10.1186/1532-429X-18-S1-O19

Cite this article as: Yang *et al.*: Towards reliable myocardial blood-oxygen-level-dependent (BOLD) CMR using late effects of regadenoson with simultaneous ¹³n-ammonia pet validation in a whole-body hybrid PET/MR system. *Journal of Cardiovascular Magnetic Resonance* 2016 **18** (Suppl 1):O19.

**Submit your next manuscript to BioMed Central
and take full advantage of:**

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at
www.biomedcentral.com/submit

