Journal of Cardiovascular Magnetic Resonance

MODERATED POSTER PRESENTATION



Detection of 3D Cardiac metabolism after injection of hyperpolarized [1-13C]pyruvate

Francesca Frijia^{1*}, Florian Wiesinger², Maria Filomena Santarelli³, Vincenzo Positano¹, Jan Henrik Ardenkjaer-Larsen⁴, Rolf Schulte², Luca Menichetti³, Giulio Giovannetti³, Giacomo Bianchi⁵, Giovanni Aquaro¹, Vincenzo Lionetti⁵, Daniele De Marchi¹, Alessandra Flori⁶, Luigi Landini¹, Fabio Anastasio Recchia⁵, Massimo Lombardi¹

From 2011 SCMR/Euro CMR Joint Scientific Sessions Nice, France. 3-6 February 2011

Introduction

MRI with hyperpolarised ¹³C represents a promising modality for *in-vivo* spectroscopy and provides a unique opportunity for *non-invasive* assessment of cardiac regional metabolism.

Purpose

We present a method based on a volumetric IDEAL spiral CSI acquisition for obtaining spatial information on the metabolism of the whole heart after intravenous injection of hyperpolarized [1-13C]pyruvate in a large animal model with a clinical 3T scanner.

Methods

Three healthy male mini-pigs (38±2 kg) were maintained under deep sedation; a dose of 20 mL of 230 mM [1-13C]pyruvate was administered over about 10 s by manual injection. Animal experiments were performed on a 3T GE Signa HDx scanner with a 13C quadrature birdcage coil.

[1-13C]pyruvate was polarized using a HyperSense DNP polariser with subsequent dissolution. The final injection solution contained 230 mM sodium [1-13C] pyruvate, 100 mM TRIS buffer, 0.27 mM Na₂EDTA and 20 μ M Dotarem with T \approx 37°C and pH \approx 7.6.

Anatomical reference images were acquired in the axial plane with standard FIESTA sequence (body coil FOV=30x30 cm², FA= 20° , TE/TR=3.8ms/7.52ms, matrix 224x160, slice thickness 5 mm, 20 slices).

Metabolic information covering the heart were obtained using a 3D IDEAL spiral CSI prescribed on the same region imaged by the reference anatomical sequence (FOV= 30x30 cm, slab thickness=100mm) starting 20 seconds after the beginning of the hyperpolarized $[1^{13}C]$ -pyruvate injection. The IDEAL spiral CSI concept was implemented into a multi-slice, pulse-and-acquire sequence with a 2D spiral readout and phase encoding along the third dimension. A constant echo time shift of TE=0.9ms, 11 encoding steps and FA=7° were used to optimize the study for the considered frequencies.

The data was reconstructed using spectrally-preconditioned, minimum-norm CS inversion followed by gridding reconstruction implemented in Matlab. The reconstruction on cardiac short axis (SA) and image fusion was performed by PMOD software.

Results

Pyruvate and its metabolic products lactate and bicarbonate were detected in the heart. Metabolic maps overlaid on anatomical images are shown in Figure 1. On SA sections the metabolites signal resulted correctly localized in cardiac structures: pyruvate more evident in ventricular cavity, bicarbonate in myocardial wall.

Conclusions

This study demonstrated that the volumetric spatial localization of ¹³C metabolites on the whole pig heart can be acquired with a 3D IDEAL spiral CSI sequence. This allows investigation of the metabolic behaviour over all the heart, both in normal and pathological conditions.

Full list of author information is available at the end of the article



© 2011 Frijia et al; licensee BioMed Central Ltd. This is an open access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/2.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

¹Fondazione G.Monasterio CNR-Regione Toscana, Pisa, Italy

Author details

¹Fondazione G.Monasterio CNR-Regione Toscana, Pisa, Italy. ²GE Global Research, Munich, Germany. ³Institute of Clinical Physiology, National Research Council, Pisa, Italy. ⁴GE Healthcare, Huginsvej, 3400 Hillerod, Denmark. ⁵Sector of Medicine, Scuola Superiore Sant'Anna, Pisa, Italy. ⁶Scuola Superiore Sant'Anna, Pisa, Italy.

Published: 2 February 2011

doi:10.1186/1532-429X-13-S1-M2

Cite this article as: Frijia et al.: Detection of 3D Cardiac metabolism after injection of hyperpolarized [1-13C]pyruvate. Journal of Cardiovascular Magnetic Resonance 2011 13(Suppl 1):M2.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

BioMed Central

Submit your manuscript at www.biomedcentral.com/submit